

Cluster size determination using shadowgraphy measurements

Hanna Eick, Philipp Brand, Christian Mannweiler, Sophia Vestrick, and Alfons Khoukaz

Institute for Nuclear Physics, WWU Münster


In cooperation with the group of O. Willi, HHU Düsseldorf

30th Conference of the International Nuclear Target Development Society INTDS, Switzerland
September 25 – 30, 2022




Motivation of Measurements

- Cluster-jet targets are important and central component of several experiments (e.g., for laser induced proton acceleration → ARCTURUS, or hadron physics → PANDA experiment)
- Investigation of target properties and cluster generation process
- Important properties of the targets are the sizes of clusters and the size distribution
 - Prediction of the beam's time structure
 - Understanding of the residual gas




Important Experiments with Cluster-Jet Targets

Laser induced proton acceleration



ARCTURUS TW laser at HHU Düsseldorf

Hadron physics




See talk by S. Vestrick

PANDA experiment at future HESR

And other experiments: See talk by P. Brand

Set-up of Münster Cluster-Jet Targets

- Hydrogen is flowing through purification cartridge and cooled down by two-stage cold head
- Insulation vacuum chamber with turbo pump to ensure a thermal decoupling
- Cooled gas is pressed through laval nozzle with diameter of $37\ \mu\text{m}$
- Depending on target conditions, different cluster production processes occur






Generated by D. Veith, edited by L. Leßmann

Cluster Production Process

- Laval nozzle with convergent inlet and divergent outlet
- Hydrogen can be either gaseous or liquid in front of the nozzle → leads to different cluster formation processes




Cluster size determination using shadowgraphy measurements



Set-up for Shadowgraphy Measurements


- Top view of the experimental set-up
- Cluster-Jet Target, developed and built up in Münster
- Target setting determines:
 - Target beam thickness
 - Cluster size distribution



Created by Christian Mannweiler

Set-up for Shadowgraphy Measurements


- ARCTURUS TW laser system of HHU Düsseldorf
- Ultrashort-pulse laser (30 fs) is used as background lightening
- Pictures of clusters are taken with a camera in combination with a microscope objective
- With longer exposure time clusters at about (200 - 1000) m/s would no longer be recognizable as dots



Created by Christian Mannweiler


Shadowgraphy Image

- Shadowgraphy image (background subtracted) with ...
 - some sharp clusters
 - a lot of interference rings
 - dust on camera and optics (some of this is eliminated by background subtraction)




Shadowgraphy Image

- Shadowgraphy image (background subtracted) with ...
 - some sharp clusters
 - a lot of interference rings
 - dust on camera and optics (some of this is eliminated by background subtraction)




Shadowgraphy Image

- Shadowgraphy image (background subtracted) with ...
 - some sharp clusters
 - a lot of interference rings
 - dust on camera and optics (some of this is eliminated by background subtraction)



Shadowgraphy Image

- Shadowgraphy image (background subtracted) with ...
 - some sharp clusters
 - a lot of interference rings
 - dust on camera and optics (some of this is eliminated by background subtraction)







Cluster Selection Criteria

- For every target setting and position in the cluster beam 1000 pictures were taken → impossible to analyze all by hand
- An automated procedure is required to find potential cluster candidates and to decide with cluster selection criteria whether they are clusters or not
- The image processing program *ImageJ*^[2] is used to identify potential cluster candidates




Cluster Selection Criteria

- Widths of both profiles must be similar in size
- Amplitude of peak must have significant height
- Background must not be too uneven
- Sorted out when light spot (Poisson spot) appears in center → too far away from focal plane




Cluster Selection Criteria

- Widths of both profiles must be similar in size
- Amplitude of peak must have significant height
- Background must not be too uneven
- Sorted out when light spot (Poisson spot) appears in center \rightarrow too far away from focal plane




Light spot




Calculation of Cluster Diameter

- After identifying real clusters, a 2-dimensional fit is applied which includes the diameter as fit parameter
- By means of a calibration measurement (wire measurement) the camera pixel information is converted into a size information in μm




Cluster Size Distribution

- Cluster size distribution (1000 pictures) for one target setting
- Steeper edge on left side and longer tail on right side → ‘real’ distribution on right side, but resolution limit on left side




Do we really see the True Distribution?



- To measure depth of every picture, calibration measurements using μm -sized toner particles were performed
- Toner inks were moved at location of target through laser in micrometer steps
- First intention: Look at one specific toner particle and find out in which range it can be found and analyzed

Do we really see the True Distribution?




- The observation made is that the ‘toner clusters’ look larger the further they are away from the focal plane
- The found distribution is not the true distribution
- Bringing together information from the measured distribution and the toner clusters (deconvolution)


True Cluster Size Distribution

Result with assumption that true distribution is symmetric Gaussian

→ More quantities can be determined (volume density, gas flow) and compared with measured values and theoretical calculations


Symmetric Gaussian







Volume Density

- For different positions in the cluster beam, a volume density (blue) can be calculated from size distributions
- Expected volume density (orange) is higher than calculated → Also smaller clusters have an influence on volume density
- Structure in cluster beam → Core beams





Outlook

- It is very likely that there are also smaller clusters, but they cannot be found with pure shadowgraphy method (at the mentioned conditions)
- Methods to find the distribution of smaller clusters will be tested in the future (3-WEM measurements)



Summary

- Shadowgraphy measurements were performed at the ARCTURUS laser in Düsseldorf to estimate the cluster size distribution of a Münster Cluster-Jet target
- A routine was developed to find, select and analyze the clusters, which can be seen as dark spots on the shadowgraphy images
- Preliminary cluster size distribution was calculated, and the information of the toner measurements is used to find out the ,real' distribution



Thanks a lot for your attention!